




 Features & Benefits



MIP recognizes that manual repositioning exposes caregivers to increased risk of neck, shoulder, and back injuries. That's why we have developed the Swift® UltraSlide System, a cost-effective, low-friction repositioning system intended to allow caregivers to safely and easily perform manual handling tasks with little physical strain. Implementing the standard use of this system will create a safer healthcare environment for caregivers, ensure the comfort and security of patients/residents, and assist in supporting your facility's no-lift policies.

- Designed and developed by MIP in partnership with a practicing Physiotherapist
- It's Ultra Thin breathable fabrics make the device suitable to remain under the patient at all times, without any negative impact on the performance of specialty mattresses
- Swift provides not only an appropriate tool for caregivers to safely reposition, boost or turn a patient but also provides comfort and dignity to the patient

The Swift® Slider is a friction reducing device that remains in-situ to assist with patient/resident repositioning, boosting and turning in bed. Made with a breathable, low-friction backing it allows staff to reposition patients/residents safely and easily whilst avoiding heat build up under the patient/resident. The Swift® Slider's thin profile will not impair the performance of pressure redistribution mattresses and can therefore remain on the bed when not in use (tucked in), making it readily available when needed.

- Standard Swift Slider - designed for both standard and bariatric use. The central panel should support patient/resident from neck to just below the buttocks
- Full Length Swift Slider - designed for both standard and bariatric use. The Full Length Swift Slider supports the full length of a patient/ resident. When the slider is tucked in (not in use) the Slider has the visual appeal of standard bedding.





The Swift UltraSlide System uses the UltraSlide® Bottom Sheet in conjunction with the Swift® Sliders for use with more dependent patients/residents who have little or no mobility and require extensive or total assistance for repositioning.

For all repositioning and manual handling solutions it is recommended to seek advice / approval from a manual handling expert.

NOTE: Inspect before use for signs of wear-and-tear and replaced as necessary. Fabric, thread and stitch degradation will take place over time (as normal wear and tear) that will reduce the slide capacity.

Specifications

STYLE	BACKING COLOUR	DESCRIPTION	SIZE
-------	----------------	-------------	------

STYLE	BACKING COLOUR	DESCRIPTION	SIZE
PTD-56/BA/Y/UK	Yellow Taffeta	Swift® Slider	123 x 186cm - Suitable for standard and bariatric application
PTD-56/BA/LG/UK	Blue	Swift Slider Long Length	L200 x W205cm - Suitable for standard and bariatric application